



System One

## BiPAP S/T



## Compact, lightweight and easy to operate

Aligned with the same platform as our System One sleep therapy technology, the Philips Respironics **BiPAP S/T noninvasive ventilator** offers a straightforward user interface and integrated heated humidification for easy setup and management. Bi-level S (spontaneous) and S/T (spontaneous/timed) triggers respond on demand or automatically to deliver prescribed therapy while System One Heated Tube humidification allows humidity levels to be tailored to meet patient needs.

### Enhanced functionality

- Increased pressure capabilities to 30 cm H<sub>2</sub>O
- Onboard memory stores patient therapy data without an SD card
- System One alignment leverages accessories across all Philips Respironics sleep therapy and noninvasive ventilator devices
- Heated Tube technology tailors humidity levels to patient needs

### System One features

- Advanced event detection
- System One Resistance Control
- Oximetry module
- Encore and DirectView patient management systems
- SD card and wired, wireless, and broadband modem capable

### Advanced leak-sensing technology with Digital Auto-Trak

- Enables optimum patient/ventilator synchrony
- Recognizes and compensates for leaks
- Automatically adjusts variable trigger and cycle thresholds

### Heated Tube humidification

- Air temperature control
- Improved humidification
- Rainout protection



As your Allies in Better Sleep and Breathing, we lead the way by providing more than advanced technology: we provide exceptional support.

**PHILIPS**  
RESPIRONICS

# BiPAP S/T ordering information

Specifications		Device ordering information	Part number
Pressure range	EPAP 4-30 cm H <sub>2</sub> O IPAP 4-30 cm H <sub>2</sub> O CPAP 4-20 cm H <sub>2</sub> O	BiPAP S/T, C Series	IN1061X
Breath rate	0-30 BPM	<b>Ventilation accessories</b>	
Flex pressure relief	0 to 3 (S-Mode)	<b>Humidification</b>	
Ramp time	0 to 45 min (5-min increments)	System One 60 Series humidifier	IN6H
Starting ramp	4 to EPAPmin	System One 60 Series Heated Tube humidifier (includes 80 watt power supply and Heated Tube)	IN6T
Dimensions	7" L x 5.5" W x 4" H 18 cm L x 14 cm W x 10 cm H	Water chamber (includes both 50 and 60 Series lids)	1063785
Weight	3 lbs; 1.36 kg (without humidifier)	Heated Tube humidifier upgrade kit	6T-CONVKIT
Filters	Foam and optional ultra-fine	<b>Tubing</b>	
Device setup	LCD/control wheel/push button	System One Heated Tube	SYS1HT15
Data storage capacity (minimum)	Display: compliance VIC, 1-, 7- and 30-day averages SD card > 1 year On-board: up to 3 months	System One Performance Tubing (15mm)	1074400
Compliance meter	Breathing detection	System One Performance Tubing (22mm)	1032907
Altitude compensation	Automatic	<b>Filters</b>	
Electrical requirements	100-240 VAC, 50/60 Hz	Pollen (foam) filter, reusable (1/pack)	1035443
Humidification	Integrates with System One 60 Series humidifier and System One 60 Series Heated Tube humidifier	Pollen (foam) filter, reusable (2/pack)	1029330
Miscellaneous	Lighted LEDs	Ultra-fine filter, disposable (1/pack)	1063096
		Ultra-fine filter, disposable (2/pack)	1035442
		Ultra-fine filter, disposable (6/pack)	1029331
		<b>Accessories</b>	
		Device side panel	1099592
		Carrying case	1063857

\*Philips Respironics devices with advanced event detection detect, respond to (some models), and log these advanced events: Apnea Hypopnea Index (AHI), respiratory-related arousal (RERA), snore, flow limitation (FL), leak, vibratory snore (VS), obstructive airway apnea (OA), large leak (LL), clear airway apnea (CA), periodic breathing (PB), and hypopnea (H).

Respironics, BiPAP, AVAPS, Bi-Flex, Encore and System One are trademarks of Koninklijke Philips N.V. and its affiliates. All rights reserved.

For more information, please visit [www.philips.com/respironics](http://www.philips.com/respironics)



© 2013 Koninklijke Philips N.V. All rights are reserved.

Philips Healthcare is part of Royal Philips

Philips Healthcare reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.

Broudy AE 7/26/13 MCI 4105741 PN 1111569

Philips Respironics Asia Pacific  
+65 6882 5282  
Philips Respironics Australia  
+61 (2) 9947 0440  
Philips Respironics Benelux  
+33 1 47 28 30 83  
Philips Respironics Deutschland  
+49 8152 93 06 0  
Philips Respironics France  
+33 2 51 89 36 00  
Philips Respironics Iberia  
+34 91 566 90 00

Philips Respironics Italy  
+39 039 203 1  
Philips Respironics Sweden  
+46 8 120 45 900  
Philips Respironics Switzerland  
+41 6 27 45 17 50  
Philips Respironics United Kingdom  
+44 800 1300 845

Philips Respironics  
1010 Murry Ridge Lane  
Murrysville, PA 15668 USA