

Dorma 200

INV201, INV201H



PHILIPS
RESPIRONICS

Specifications

Pressure Range:	4 to 20 cm H ₂ O
Flex pressure relief:	On or off
Ramp time:	0 to 45 min (5-min increments)
Starting ramp pressure	Patient adjustable
Dimensions:	7" L x 5.5" W x 4"H 18 cm L x 14 cm W x 10 cm H
Weight:	3lbs; 1.36 kg (without humidifier)
Filters:	Foam and optional ultra-fine
Device set-up	LED/keypad
Data Storage Capacity:	Compliance VIC 7-30-day averages On-board: > 3 years SD Card (optional) : Lifetime
Compliance meter:	Breathing detection
Altitude compensation:	Automatic
Electrical Requirements:	100-240 VAC, 50/60 Hz
Humidification :	Integrates with System One heated humidifier
Humidification Modes:	Classic
Warranty:	Two years
Miscellaneous:	Lighted LEDs