

Pulmo-Aide® Compressor Nebulizer System

The Gold Standard for High Performance Medication Delivery



Pulmo-Aide Compressor
Nebulizer System
Model 5650D

The DeVilbiss Pulmo-Aide Compressor Nebulizer System has been providing **effective aerosol treatments for more than thirty years**. Modern styling, durability and ease of use make this unit **perfect for in-home asthma treatments**. The Pulmo-Aide System comes ready to use with all the attachments – filter, tubing, and nebulizer.

The **Pulmo-Aide System is the only approved compressor for use with TOBI** (tobramycin solution for inhalation). The Pulmo-Aide is unmatched in quality, value and performance.

Performance Features

- Maximum compressor pressure – 30 psig
- Compressor free air flow – 12 lpm
- 5-year limited warranty
- Longest lasting compressor on the market with an average product life of 10 years
- Unmatched quality and durability – Great for nursing homes and long-term users

Design Features

- Metal frame motor construction
- Unique 4-pole motor design extends life of the unit
- Large built-in compartment for storing nebulizers and tubing
- Built-in nebulizer holder is a convenient place to store your nebulizer if your treatment is interrupted.

Available Models

5650D	Pulmo-Aide Compressor Nebulizer System 115VAC, 2 prong plug
5650C	Pulmo-Aide Compressor Nebulizer System CSA approved

Parts and Accessories

DEVILBISS PART #	PRODUCT DESCRIPTION
900D-610	VixOne Reusable Nebulizer (case of 10)
3655D-621	VixOne Disposable Nebulizer (case of 50)
3655PN-612	PediNeb™ Pacifier (12 pack)
AAM12	Adult Aerosol Mask (12 pack)
AAM20	Adult Aerosol Mask (50 pack)
PAM10	Pediatric Aerosol Mask (50 pack)
DL1050	Pediatric Character Mask (50 pack)
5650D-601	Compressor Filter (package of 25)
5650D-602	Compressor Filter (package of 5)
5650D-612	Filter Kit (6 each – caps and filters)
5650D-611	Tubing Connector (package of 6)

VixOne™ Nebulizer

- MMAD of 2.7 microns
- Respirable Fraction – 84% below 5 microns
- 7' kink-resistant flexhose tubing
- Latex-free
- Reusable model is dishwasher safe to promote frequent cleaning
- Low residual volume (0.7cc)
- High capacity medication bowl (>10cc)
- Comfortable mouthpiece design
- Anti-spill free design provides secure seal
- Nebulizes in upright position or at angles up to 45° to enhance comfort during treatment
- Universal connector accommodates standard adult and pediatric masks, as well as the PediNeb pediatric pacifier mouthpiece

Technical Specifications

Built to meet the requirements of HCPCS code: E0570		Warranty:	5 Years
Product Dimensions:	10.5"H X 6.3"W X 10.1"D	Safety Approvals:	ETL
Product Weight:	7.1 lbs.	Motor/Pump Assembly:	Diaphragm, Primarily Metal Parts
Maximum Air Pressure:	30 psig*	Neb Storage:	Yes
Maximum Liter Flow:	12 lpm*	Master Carton Quantity:	4
MMAD:	5 microns or less		
Sound:	51-55 dBA*		
Filter Replacement:	No tools required		

*Specifications are typical and may vary slightly from system to system



DeVilbiss Healthcare LLC | 100 DeVilbiss Drive | Somerset, PA 15501 USA
800.338.1988 | 814.443.4881 | www.DeVilbissHealthcare.com