



Kangaroo™ 924
Enteral Feeding Pump





When ease of use, accuracy and dependability matter, the KANGAROO™ 924 Enteral Feeding Pump delivers it all:

- Touch panel and knob operation for easy-to-use set up and programming.
- Preset volume and over-infusion safeguard ensure correct and consistent delivery of nutrition.
- Audible and visual alarms reinforce clinician confidence and patient safety.
- Highly visible LED display and built-in pole clamp aid in ease of use.

Features

- Accurate pre-set volumetric delivery with volume totaliser
- Intuitive and user-friendly controls
- Easy-to-read LED display
- Built-in pole clamp
- Over-infusion safeguard
- Audio and visual alarms
- 24-hour memory
- 1-year warranty

Specifications

- Pump Mechanism: Rotary Peristaltic Mechanism
- Pump Set Required: KANGAROO™ 924 Pump Sets (refer below)
- Operating Range: 1 - 300 ml/hour
- Volume to be Delivered: 1-2000 ml in 1ml/hour increments
- Accuracy: $\pm 10\%$
- Occlusion Pressure: 12 psi
- Alarms: Flow Error, Low Battery, Hold Error, No Set, System Error, Volume To Be Delivered (VTBD)
- Battery: Nickel Metal Hydride, Rechargeable
- Battery Life: 16 hours
- Memory: 24 hours
- Dimensions: 18.4 cm (High) x 14 cm (Wide) x 10.8 cm (Deep)
- Weight: 1.64 kg
- Languages: English, Spanish, French, Portuguese, Dutch, Danish, German, Italian, Norwegian, Swedish, Turkish, Finnish, Polish



Proximal Spike



3-In-1 Set



1000mL Bag

Code	Description	Qty
KANGAROO™ 924 Enteral Feeding Pump		
492457	924 Enteral Feeding Pump	1
Pump Set (Sterile Sets)		
777100	1000 mL, Sterile, Bag	30
777050	500 mL, Sterile, Bag	30
777301	3-in-1 Spike Set, Sterile	30
674626	Proximal Spike Set	30
673622	1000 mL Pump Set, DEHP-Free	30



Kangaroo™ 924 Pump Set



COVIDIEN

positive results for life™

COVIDIEN, COVIDIEN with Logo "positive results for life" and ™
marked brands are trademarks of Covidien AG or its affiliate.
© 2009 Covidien AG or its affiliate. All rights reserved.

Ent 037-02-09

TYCO HEALTHCARE PTY LTD
166 EPPING ROAD,
LANE COVE NSW 2066
AUSTRALIA
(T) 1800 252 467

TYCO HEALTHCARE LTD
GROUND FLOOR, 15B VESTEY DRIVE,
MOUNT WELLINGTON, AUCKLAND
NEW ZEALAND
(T) 0508 489 264

WWW.COVIDIEN.COM