



Quality for a  
reasonable price

# PARI COMPACT

- Inhalation device for diseases of the lower airways, such as **COPD, Asthma or Bronchitis**
- **Proven efficacy:** includes PARI LC PLUS® nebuliser with mouthpiece and child's mask

## Technical data:

Operating pressure: 1.0 bar

Weight: 1.2 kg

Dimensions (w x h x d) 10 x 14 x 16 [cm]

Warranty: 2 Years

## Aerosol Performance:

Total Output Rate: 418 mg/min

MMD: 3.9 µm (Mass Median Diameter)

Mass percentage below 5 µm: 64 %

Measured with the Malvern MaserSizer X (according to Fraunhofer-Model) at 21°C and 50% relative humidity.  
Nebulised medium: 0.9% NaCl (5 ml). Inspiratory flow 20 l/min

Specialists in effective inhalation



# PARI COMPACT

and its accessories



## PARI COMPACT

Inhalation device  
230-240 V~/50Hz  
With PARI LC PLUS® Nebuliser,  
Child Mask, Mouthpiece uni-  
versal and Connecting tubing

Item No. 052G1002



## Connection Tubing (f/f)

Item No. 041B4592



## PARI LC PLUS® Year Pack

consists of LC PLUS® Nebuliser,  
Connection Tubing and Filter

Item No. 022G8010



## PARI Filter Holder with Filter

Item No. 041B6500



## PARI Child Mask Soft „Spiggy“

with Elastic Band

Item No. 041G0741



## PARI Adult Mask Soft

with Elastic Band

Item No. 041G0470



PARI COMPACT  
Item No. 052G1002

**PARI**  
Export Service Center  
**+49 (0) 8151 279-220**

Give us a ring ...

**PARI** GmbH  
Specialists in effective inhalation  
Moosstrasse 3 · D-82319 Starnberg  
Telefon +49(0) 8151 279-0  
Fax +49(0) 8151 279-101

You can find PARI on the internet:  
E-Mail: [info@pari.de](mailto:info@pari.de) · [www.pari.de](http://www.pari.de)

„misprint excepted“

Specialists in effective inhalation

